

Semi Automated Plate Sealer

Features:

- * Format different micro well plates and heat seals
- * Adjustable Sealing Temperature: 80 - 200°C
- * OLED display screen, high light and no visual angle limit
- * Plate adapters allow use of PCR plate or virtually any ANSI micro plate.
- * Sealing time and temperature can be set by using the control knob, OLED display screen, high light and no visual angle limit.


Control:

The MPSE22 Semi Automated Plate Sealer allows optimal settings for all types of plates and sealing material. Sealing time and temperature can be set using the control knob.




- * Auto
- * Sealing Time in 0.1 sec increments
- * Sealing Temperature in 1.0°C increments

Accessories:



96-well PCR plate adapter



Standard plate adapter



Weighted sealing flake

Technical Parameters:

Model	MPSE22
Adjustable Sealing Temperature	80~200°C
Temperature Accurate	1.0°C
Seal Time	0.5s~10s, increased by 0.1 sec.
Seals Plate Heights	9 to 48mm
Compatible Plate Materials	PP(Polypropylene),PS(Polystyrene),PE(Polyethylene)
Compatible Plate Types	SBS standard plates, Deep-well plates, PCR plates(Skirted, semi-skirted and no-skirted formats)
Consumption	300W
Power Supply	AC110/220V±10%,50/60Hz
External Size (L*D*H)	370*178*330mm
Package Size(L*D*H)	530*320*540mm
Net/Gross Weight	9.6/12kg