

## Vacuum Pump

**VanPel**  
Laboratoriumapparatuur



### Features:

- \* Oil free, clean exhausted air by built-in filter.
- \* Automatic cooling exhaust system ensures 24h operation continuously.
- \* The pressure adjustable design can satisfy a certain range of vacuum and gas velocity.
- \* TEFLON surface treatment, anti-gas corrosion.
- \* Orange look for anticorrosive vacuum pump.(eg: GM-0.33IIP, GM-0.5P, GM-0.5IIP, GM-1.0P)

### Technical Parameters:

Model	VAC-12	VAC-30-1	VAC-30-2	VAC-20	VAC-60
Pump Head	1	2	1	1	2
Speed of Evacuation	12L/Min	30L/Min	30L/Min	20L/Min	60L/Min
Ultimate Pressure	≥0.075Mpa	≥0.095Mpa	≥0.08Mpa	≥0.08Mpa	≥0.08Mpa
Vacuum	250mbar	50mbar	200mbar	200mbar	200mbar
Pressure	≥30 Psi Positive or negative pressure	Negative pressure	≥30 Positive, Positive or negative pressure	Negative pressure	≥30 Positive, Positive or negative pressure
Inlet	φ6mm	φ6mm	φ6mm	φ6mm	φ6mm
Outlet	φ6mm	Silencer	φ6mm	Silencer	φ6mm
Working Temp.	7-40°C				
Body Temp.	<55°C				
Noise	<60dB			<50dB	
Motor Power Consumption	75W	160W			
Power Supply	AC110/220V±10%, 50/60Hz				
Diaphragm	NBR				
Valves	/	NBR			
Net Weight(kg)	4	10	6.6	6.2	10
Gross Weight(kg)	4.5	11	7.6	7.5	11
External Size(L*W*H) mm	195*98*156	300*120*235	210*160*235	215*120*235	300*160*235
Package Size(L*W*H) mm	330*210*240	314*205*316	300*210*310	315*200*310	380*340*270