

## Vacuum Drying Oven



### Features:

- \* PID microprocessor temperature control with LED display.
- \* Polished stainless steel inner chamber.
- \* Vacuum degree is automatically controlled
- \* Silicon rubber seal ring ensures the high vacuum degree in the chamber.
- \* Tempered double-layer glass door ensures the operation safety.
- \* The vacuum is controlled by microprocessor controller with digital display
- \* The working program can be cycled for 99 times
- \* The vacuum pump is standard accessory
- \* Inert gas can be filled in the chamber. (Inert gas pressure is not higher than 0.1Mpa when filling in)

### Technical Parameters:

Model	VAC30M	VAC50M	VAC90A	VAC215A	VAC90AD	VAC215AD
Capacity	30L	54L	91L	215L	91L	215L
Temp. Range	50~200°C	50~250°C	50~200°C			
Temp. Fluctuation	±1°C					
Vacuum Degree	<133Pa					
Vacuum Control	/			Automatically; Include one vacuum meter		Automatically; Microprocessor controller with digital display
Shelves No.	1 pc	2 pcs	2 pcs	3 pcs	2 pcs	3 pcs
Consumption	500W	1400W	1600W	2200W	1600W	2200W
Power Supply	AC110V/220V±10%, 50/60Hz					
Internal Size(W*D*H) mm	330*330*280	420*350*370	450*450*450	560*600*640	450*450*450	560*600*640
External Size(W*D*H) mm	580*570*450	710*590*530	610*660*1450	720*820*1750	610*660*1450	720*820*1750
Package Size(W*D*H) mm	700*640*600	820*700*690	790*750*1570	900*850*1930	790*750*1570	900*850*1930
Gross Weight(kg)	50	106	155	195	155	195